

SureStart *Taq* DNA Polymerase

INSTRUCTION MANUAL

Catalog #600280 (100 U), #600282 (500 U), #600284 (1000 U)

Revision A.01

For in Vitro Use Only

600280-12

LIMITED PRODUCT WARRANTY

This warranty limits our liability to replacement of this product. No other warranties of any kind, express or implied, including without limitation, implied warranties of merchantability or fitness for a particular purpose, are provided by Agilent. Agilent shall have no liability for any direct, indirect, consequential, or incidental damages arising out of the use, the results of use, or the inability to use this product.

ORDERING INFORMATION AND TECHNICAL SERVICES

United States and Canada

Agilent Technologies

Stratagene Products Division

11011 North Torrey Pines Road

La Jolla, CA 92037

Telephone (858) 373-6300

Order Toll Free (800) 424-5444

Technical Services (800) 894-1304

Internet techservices@agilent.com

World Wide Web www.stratagene.com

Europe

Location	Telephone	Fax	Technical Services
Austria	0800 292 499	0800 292 496	0800 292 498
Belgium	00800 7000 7000	00800 7001 7001	00800 7400 7400
	0800 15775	0800 15740	0800 15720
France	00800 7000 7000	00800 7001 7001	00800 7400 7400
	0800 919 288	0800 919 287	0800 919 289
Germany	00800 7000 7000	00800 7001 7001	00800 7400 7400
	0800 182 8232	0800 182 8231	0800 182 8234
Netherlands	00800 7000 7000	00800 7001 7001	00800 7400 7400
	0800 023 0446	+31 (0)20 312 5700	0800 023 0448
Switzerland	00800 7000 7000	00800 7001 7001	00800 7400 7400
	0800 563 080	0800 563 082	0800 563 081
United Kingdom	00800 7000 7000	00800 7001 7001	00800 7400 7400
	0800 917 3282	0800 917 3283	0800 917 3281

All Other Countries

Please contact your local distributor. A complete list of distributors is available at www.stratagene.com.

SureStart *Taq* DNA Polymerase

CONTENTS

Materials Provided	1
Storage Conditions	1
Notice To Purchaser	1
Introduction	2
Critical Optimization Parameters	3
Primer–Template Purity and Concentration	3
Enzyme Concentration	3
Cycling Parameters.....	3
Activating SureStart <i>Taq</i>	4
Primer Design.....	5
Minimizing the Effects of Contaminating DNA	5
Reaction Buffer	5
Polymerase Chain Reaction Protocol Using SureStart <i>Taq</i> DNA Polymerase	6
Setting Up the PCR Reaction	6
Examination of the PCR Products Using Gel Electrophoresis.....	7
Troubleshooting	8
References	9
MSDS Information	9

SureStart Taq DNA Polymerase

MATERIALS PROVIDED

Material Provided	Quantity		
	#600280	#600282	#600284
SureStart Taq DNA polymerase (5 U/ μ l)	100 U	500 U	1000 U
SureStart Taq 10 \times reaction buffer	1 ml	2 \times 1 ml	4 \times 1 ml

STORAGE CONDITIONS

All Components: -20°C

NOTICE TO PURCHASER

Notice to Purchaser: Limited License

Purchase of this product includes an immunity from suit under patents specified in the product insert to use only the amount purchased for the purchaser's own internal research. No other patent rights (such as 5' Nuclease Process patent rights) are conveyed expressly, by implication, or by estoppel. Further information on purchasing licenses may be obtained by contacting the Director of Licensing, Applied Biosystems, 850 Lincoln Centre Drive, Foster City, California 94404, USA.

Revision A.01

© Agilent Technologies, Inc. 2009

INTRODUCTION

Stratagene SureStart *Taq* DNA polymerase is a modified version of *Taq2000* DNA polymerase with hot start capability. SureStart *Taq* DNA polymerase improves PCR amplification reactions by decreasing background and increasing amplification of desired products. Using SureStart *Taq*, hot start is easily incorporated into PCR protocols already optimized with *Taq* DNA polymerase, with little or no modification of cycling parameters or reaction conditions.

Preparing PCR reaction mixtures at room temperature can result in high background and lower yields of specific products. Certain PCR enzymes exhibit significant polymerase activity at temperatures encountered during reaction setup or while ramping up to stringent primer annealing temperatures. For example, *Taq* DNA polymerase exhibits 2–3% maximum activity at 25°C (e.g., room temperature setup) and 70% maximum activity at 50°C (which is generally below the melting temperature of PCR primers).¹ Nonspecific primer annealing and extension at nonrestrictive temperatures produces undesirable products that are amplified throughout the remaining cycles. Misprimed products and artifacts such as primer-dimers can impair gel analysis, quantitation, and sequencing of specific PCR products. In the amplification of misprimed products, dNTPs and primers are diverted away from the synthesis of specific products, reducing overall yields.

A number of hot start techniques have been developed to improve amplification specificity and to allow PCR setup at ambient temperatures.^{2–6} The most convenient hot start methods employ reversibly inactivated enzymes. With SureStart *Taq* DNA polymerase, hot start is provided by a modified form of Stratagene *Taq2000* DNA polymerase, a highly purified recombinant version of *Taq* DNA polymerase. SureStart *Taq* DNA polymerase remains inactive until stringent temperatures (e.g., 92–95°C) are reached. SureStart *Taq* can be activated by adding an initial step of 9–12 minutes at 92–95°C to the beginning of PCR cycling programs. Alternatively, the enzyme can be activated slowly during temperature cycling without prior activation, although it may be necessary to add additional cycles to existing cycling programs to achieve optimal product yield. Either activation method provides a PCR hot start, since primer extension can not occur during PCR setup when SureStart *Taq* DNA polymerase is inactive. SureStart *Taq* can be used in a variety of amplification systems, including quantitative PCR and RT-PCR, to improve specificity, yield, and amplification of difficult targets.

SureStart *Taq* DNA polymerase is optimized for use in amplifying DNA targets between 100 and 2000 bases. Although templates of up to 5 kb may be amplified with SureStart *Taq* DNA polymerase, we recommend the use of Stratagene Herculase enhanced DNA polymerase to achieve optimal yields of longer targets (>2 kb).

CRITICAL OPTIMIZATION PARAMETERS

SureStart *Taq* DNA polymerase can be introduced into all existing amplification systems optimized with *Taq* or *Taq2000* DNA polymerase, with only slight modification of cycling conditions. For other systems, optimization may be required to achieve highest product yield and specificity. Critical optimization parameters include the quality and concentration of DNA template and primers; the concentration of enzyme; cycling parameters; primer design; and the use of recommended reaction buffers.

Primer–Template Purity and Concentration

The most successful PCR results are achieved when amplification reactions are performed using purified primers and templates that are essentially free of extraneous salts. Gel-purified primers, generally >18 nucleotides in length, are strongly recommended for use in SureStart *Taq*-based PCR.

Additionally, an adequate concentration of primers and template should be used to ensure a good yield of the desired PCR products. When DNA of known concentration is available, the reaction should contain at least 10^4 DNA template molecules, but no more than 1 μg total template DNA. Amounts of 25–100 ng of DNA template/100- μl reaction are typically used for amplifying single-copy chromosomal targets. Amplifying a single-copy target from complex genomic DNA is generally more difficult than amplifying a fragment from a plasmid or phage. Less DNA template can be used for amplifying lambda or plasmid PCR targets or for amplifying multicopy chromosomal genes (typically 10–100 ng).⁷ The mutation frequency can be reduced by limiting the number of PCR cycles; however, a corresponding increase in DNA template concentration is required to achieve comparable yields of PCR product. We recommend using primers at a final concentration of 0.1–0.5 μM , which is equivalent to ~100–250 ng of an 18- to 25-mer oligonucleotide primer in a 100- μl reaction volume.

Enzyme Concentration

The concentration SureStart *Taq* DNA polymerase required for optimal PCR product yield and specificity depends on the individual target system to be amplified. Successful amplification can usually be achieved using 2.5 U of enzyme/100- μl reaction for PCR targets that are <2 kb. Further optimization will be required for targets greater than 2 kb.

Cycling Parameters

As with all PCR reactions, cycling parameters are critical for successful amplification and may require further optimization.

Denaturation Conditions

High denaturation temperatures can potentially damage DNA templates, so the denaturation temperature should be as low as possible. Denaturation temperatures of 92–95°C work well for most targets.

Extension Time and Temperature

Use an extension temperature of 72°C and extension times of 0.5–1.0 minute for targets ≤ 1 kb and 1 min per kb for targets > 1 kb in length. Longer extension times may produce higher yields for difficult targets, such as high-complexity targets or targets > 2 kb in length.

Activating SureStart *Taq*

Activation of the SureStart *Taq* DNA polymerase can be achieved by two methods: using a pre-PCR heat-activation method or by activating slowly over the course of the PCR amplification reaction (PCR heat-activation method).

Using the Pre-PCR Heat-Activation Method

In this method, SureStart *Taq* DNA polymerase is partially or almost completely activated prior to temperature cycling. Activation is accomplished by adding an incubation step of 9–12 minutes at 92–95°C to the beginning of cycling programs. For certain amplification systems, the use of longer initial denaturation times (11–18 minutes) may improve product yield; however, background may also increase.

Using the PCR Heat-Activation Method

For this method, the pre-PCR activation step is omitted and the enzyme is allowed to slowly activate during temperature cycling. For slow activation, we recommend using a standard initial denaturation step, e.g., 1 minute at 92–95°C. Due to lower initial DNA polymerase concentrations, product yield may be lower for some amplification systems after 30 cycles. In these cases, the use of additional cycles (5–10 cycles) is recommended to achieve optimal product yield.

The cycling conditions in Table II are recommended for use with the Pre-PCR Heat-Activation Method. When slow activation of SureStart *Taq* DNA polymerase is desired, the protocol can be modified as described in the footnotes of Table II.

Primer Design

Primer pairs that exhibit similar melting temperatures and complete complementarity to the template are recommended. Primer sequences should be analyzed for potential duplex and hairpin formation as well as false priming sites in order to obtain the highest yield of specific PCR products. Depending on the primer design and the desired specificity of the PCR amplification reaction, melting temperatures between 55° and 80°C generally yield the best results.⁷ The following formula⁸ is commonly used for estimating the melting temperature (T_m) of the primers:

$$T_m (^{\circ}\text{C}) \cong 2(N_A + N_T) + 4(N_G + N_C)$$

where N equals the number of primer adenine (A), thymidine (T), guanine (G), or cytosine (C) bases. Several other articles present additional equations for estimating the melting temperature of the primers.^{9,10} Finally, care must be taken when using degenerate primers. Degenerate primers should be designed with the least degeneracy at the 3' end. Optimization of degenerate primer concentration is necessary.

Minimizing the Effects of Contaminating DNA

To minimize the possibility of DNA contamination, use dedicated pipettors with filter-tip pipet tips. Assemble PCR reactions in an area in the lab distinct from areas where PCR analysis is taking place. Dispensing buffer, sterile water, MgCl₂, primers, and dNTPs into small aliquots prior to use is prudent as small aliquots can be thrown away if they become contaminated.

Reaction Buffer

For optimal yield and specificity it is essential that the recommended PCR buffer is used. Suboptimal results will be achieved using other buffers including the *Taq2000* DNA polymerase PCR buffer.

POLYMERASE CHAIN REACTION PROTOCOL USING SURESTART TAQ DNA POLYMERASE

Setting Up the PCR Reaction

1. Prepare a reaction mixture for the appropriate number of samples to be amplified. Add the components in order and mix gently. Table I provides an example of a reaction mixture for the amplification of a typical genomic or low complexity target. The recipe in Table I is for one reaction and must be adjusted for multiple samples. The final volume of each sample reaction is 100 μ l. Prior to setting up the reaction mixture, review *Critical Optimization Parameters* for guidelines.

TABLE I

SureStart Taq DNA Polymerase Example Reaction Mixture

Component	Volume
dH ₂ O, sterile	To 100 μ l
SureStart Taq 10 \times reaction buffer	10 μ l
dNTP mix (25mM each dATP, dCTP, dGTP, dTTP)	0.8 μ l
Template DNA	
Genomic (100 ng/ μ l) OR Low-complexity template (λ DNA or cloned DNA) (10 ng/ μ l)	1 μ l
Primer #1 (100 ng/ μ l)	2.5 μ l
Primer #2 (100 ng/ μ l)	2.5 μ l
SureStart Taq DNA polymerase (5 U/ μ l)	0.5 μ l

2. Immediately before thermal cycling, aliquot 100 μ l of the reaction mixture into the appropriate number of sterile thin-wall PCR tubes (recommended) or standard 0.5-ml microcentrifuge tubes.
3. If the temperature-cycler is not equipped with a heated cover, overlay each reaction with ~50 μ l of DNase-, RNase-, and protease-free mineral oil (available from Sigma-Aldrich, Inc.).
4. Cycle the reactions according to the guidelines in Table II.

TABLE II**Recommended Cycling Parameters^{a,b}**

Segment	Number of cycles	Temperature	Duration
1	1	95°C	9–12 minutes ^c
2	30 ^d	95°C	0.5–1 minute
		Primer $T_m - 5^\circ\text{C}$ ^e	0.5–1 minute
		72°C	0.5–1 minute for targets ≤ 1 kb OR 1 minute /kb for targets > 1 kb
3	1	72°C	10 minutes

^a The recommended cycling conditions are for reactions carried out in thin-walled tubes.

^b Cycling conditions may require modification depending on the thermal cycler used and the particular target being amplified. See cycling protocol recommendations provided by the thermal cycler manufacturer.

^c For the pre-PCR heat-activation method, incubation for 9–12 minutes at the desired denaturation temperature (92–95°C) provides optimal results for most amplification systems. For the PCR heat-activation method, use an initial denaturation step of 1 minute at 92–95°C.

^d For the PCR heat-activation method, cycle number can be increased to 35–40 cycles to achieve optimal product yield.

^e The annealing temperature may be lowered or raised further if necessary to obtain optimal results. Typical annealing temperatures will range between 55 and 60°C.

Examination of the PCR Products Using Gel Electrophoresis

Analyze the PCR amplification products by electrophoresis using an appropriate percentage acrylamide or agarose gel.

TROUBLESHOOTING

Observation	Suggestion
No PCR product or lower yield than expected	Ensure that template DNA is of sufficient quantity and quality (not sheared or nicked); increase template amount.
	Optimize the SureStart <i>Taq</i> DNA polymerase concentration; test increases in polymerase concentration of 0.5 units per 100 μ l reaction
	Optimize the MgCl ₂ concentration by increasing in increments of 0.2 mM
	For the pre-PCR heat-activation method, optimize the temperature or time used for the initial denaturation/activation step by increasing increments of 1°C or 1 minute, respectively
	For the PCR heat-activation method, increase the total number of cycles to 35-45 cycles. For the pre-PCR heat-activation method, optimize cycle number by increasing in 3 cycle increments
	Optimize annealing temperature by lowering in 2°C increments
	Optimize annealing and extension times by increasing in 30 second increments. Allow at least 30 seconds of extension time for targets up to 1 kb, and increase extension times to 1 min per kb for targets > 1 kb
	Ensure that the recommended 10 \times buffer supplied with SureStart <i>Taq</i> DNA polymerase was used
	Confirm that primers are annealing or select alternate primer pairs for testing
	Use primer pairs with matched primer melting temps (T _m) and complete complementarity to the DNA template
	Analyze the primers sequences to ensure that duplexes or hairpins do not form
	Gel-purified or HPLC purified primers \geq 18 nucleotides in length are desired for successful PCR
	Adjust ratio of primer vs. template to optimize yield
Use thin-walled tubes	
Product contains multiple bands	Optimize annealing temperatures by raising in 2°C increments
	Use primer pairs with complete complementarity to target
	Verify that the primers hybridize only to the desired sequences on the template
	Analyze the primers to ensure that duplexes or hairpins do not form
	Gel-purified or HPLC purified primers \geq 18 nucleotides in length are desired for successful PCR
	Select alternate primer pairs for testing
Try Stratagene Perfect Match PCR enhancer	
Product appears as a smear following agarose gel electrophoresis	Use clean pipettors, filter-tip pipet tips, and fresh solutions when setting up PCR reactions; work in an area in the lab distinct from areas where PCR analysis is taking place
	Optimize SureStart <i>Taq</i> DNA polymerase concentration by decreasing in increments of 0.5 units per 100 μ l reaction
	Optimize denaturation time by increasing duration in increments of 30 seconds
	Optimize denaturation temperature by increasing in 1°C increments
	Optimize extension time by reducing in 30 second increments
	Optimize cycle number by reducing in 3 cycle increments

REFERENCES

1. Nielson, K., Cline, J. and Hogrefe, H. H. (1997) *Strategies* 10(2):40-43.
2. Mullis, K. B. (1991) *PCR Methods Appl* 1(1):1-4.
3. D'Aquila, R. T., Bechtel, L. J., Videler, J. A., Eron, J. J., Gorczyca, P. *et al.* (1991) *Nucleic Acids Res* 19(13):3749.
4. Chou, Q., Russell, M., Birch, D. E., Raymond, J. and Bloch, W. (1992) *Nucleic Acids Res* 20(7):1717-23.
5. Erlich, H. A., Gelfand, D. and Sninsky, J. J. (1991) *Science* 252(5013):1643-51.
6. Kellogg, D. E., Rybalkin, I., Chen, S., Mukhamedova, N., Vlasik, T. *et al.* (1994) *Biotechniques* 16(6):1134-7.
7. Innis, M. A., Gelfand, D. H., Sninsky, J. J. and White, T. J. (1990). *PCR Protocols: A Guide to Methods and Applications*. Academic Press, New York.
8. Thein, S. L. and Wallace, R. B. (1986). *Human Genetic Diseases: A Practical Approach*. IRL Press, Herndon, Virginia.
9. Rychlik, W., Spencer, W. J. and Rhoads, R. E. (1990) *Nucleic Acids Res* 18(21):6409-12.
10. Wu, D. Y., Ugozzoli, L., Pal, B. K., Qian, J. and Wallace, R. B. (1991) *DNA Cell Biol* 10(3):233-8.

MSDS INFORMATION

The Material Safety Data Sheet (MSDS) information for Stratagene products is provided on the web at <http://www.stratagene.com/MSDS/>. Simply enter the catalog number to retrieve any associated MSDS's in a print-ready format. MSDS documents are not included with product shipments.