

MessageMap[®] Northern Blot

INSTRUCTION MANUAL

BN #775900-12

Revision #023002a

For In Vitro Use Only



LIMITED PRODUCT WARRANTY

This warranty limits our liability to replacement of this product. No other warranties of any kind, express or implied, including without limitation, implied warranties of merchantability or fitness for a particular purpose, are provided by Stratagene. Stratagene shall have no liability for any direct, indirect, consequential, or incidental damages arising out of the use, the results of use, or the inability to use this product.

ORDERING INFORMATION AND TECHNICAL SERVICES

United States and Canada

Stratagene

11011 North Torrey Pines Road

La Jolla, CA 92037

Telephone (858) 535-5400

Order Toll Free (800) 424-5444

Technical Services (800) 894-1304

Internet tech_services@stratagene.com

World Wide Web www.stratagene.com

Stratagene European Contacts

Location	Telephone	Fax	Technical Services
Austria	0800 292 499	0800 292 496	0800 292 498
			+43 (0)1795 67036
Belgium	00800 7000 7000	00800 7001 7001	00800 7400 7400
	0800 15775	0800 15740	0800 15720
France	00800 7000 7000	00800 7001 7001	00800 7400 7400
	0800 919 288	0800 919 287	0800 919 289
			+33 (0)141 91 86 42
Germany	00800 7000 7000	00800 7001 7001	00800 7400 7400
	0800 182 8232	0800 182 8231	0800 182 8234
			+49 (0)699 509 6197
Netherlands	00800 7000 7000	00800 7001 7001	00800 7400 7400
	0800 023 0446	+31 (0)20 312 5700	0800 023 0448
	+31 (0)20 312 5600		+31 (0)20 312 5792
Switzerland	00800 7000 7000	00800 7001 7001	00800 7400 7400
	0800 563 080	0800 563 082	0800 563 081
United Kingdom	00800 7000 7000	00800 7001 7001	00800 7400 7400
	0800 917 3282	0800 917 3283	0800 917 3281
			+44 (0)207 365 1056

All Other Countries

Please contact your local distributor. A complete list of distributors is available at www.stratagene.com.

MessageMap® Northern Blot

Materials Provided	1
Storage Conditions	1
Additional Materials required	1
Introduction	2
Membrane and Marker Orientation	2
Preprotocol Considerations	3
β -actin cDNA control	3
Hybridization.....	3
Stripping the Membrane	3
Protocol	4
Preparation of Labeled Probe	4
Hybridization of Radiolabeled Probe to the MessageMap® Blot	4
Hybridization Conditions	5
Washing the Blot	5
Detection	6
Stripping the Blot for Reuse	6
β -actin cDNA Control	6
Appendix I: Hybridization and Melting Temperatures	7
Hybridization Temperature	7
Calculation of the Melting Temperature	7
Troubleshooting	7
Preparation of media and reagents	8
References	9
Endnotes	9

MessageMap[®] Northern Blot

MATERIALS PROVIDED

Material Provided	Description
MessageMap [®] northern blot	6 × 8 cm membrane
MiracleHyb [™] hybridization solution	25 ml
MiracleHyb [™] probe preparation buffer	1 ml
β-actin cDNA control (100 ng purified)	25 ng/μl

STORAGE CONDITIONS

All Components: +4°C

ADDITIONAL MATERIALS REQUIRED

Prime-It[®] II Random Primer Labeling Kit (Stratagene Catalog #300385)

[α-³²P]dCTP (3000 Ci/mmol) or [α-³²P]dATP (3000 Ci/mmol)

NucTrap[®] Probe Purification Columns (Stratagene Catalog #400701 and #400702)

Revision #023002a

Copyright © 2003 by Stratagene.

INTRODUCTION

Northern blot analysis of RNA remains the preferred method for assessing the size, number, and relative abundance of gene transcripts. Each lane of the MessageMap® northern blot contains 2 µg of twice selected, high quality poly-A(+) RNA from murine or human tissues and cell lines. The poly-A(+) RNA is size fractionated on a 1% formaldehyde agarose gel and is then transferred to positively charged nylon membranes by positive-pressure alkaline transfer. MessageMap northern blots exhibit high sensitivity and low background signal making them ideal for detection of low abundance messages. MessageMap northern blots are effective when used with traditional radiolabeled probes or with nonradioactive/chemiluminescent protocols. In addition, MessageMap northern blots can be stripped and reprobbed for ease of use and increased cost effectiveness.

Membrane and Marker Orientation

Each MessageMap northern blot is configured for clear analysis of signal. The Stratagene logo is located in the lower left hand corner of each membrane. This is the 'front' of the membrane. The RNA molecular weight marker is depicted by a series of dots in the first lane on the left (see Figure 1). Marker sizes are as follows: 9.49 kb, 7.46 kb, 4.4 kb, 2.37 kb. For individual MessageMap blot configurations please refer to the Certificate of Analysis included with your order.




FIGURE 1 Schematic diagram of the MessageMap northern blot

PREPROTOCOL CONSIDERATIONS

β -actin cDNA control

The β -actin control template is provided as a control for hybridization and may also be used for relative quantitation of a mRNA of interest. Due to the abundance of β -actin in most RNA samples, Stratagene recommends probing with the β -actin cDNA control only **after** detection of the mRNA of interest. This will ensure that signal from the β -actin cDNA control does not overwhelm signal from lower abundance sample RNAs. The human cytoplasmic β -actin cDNA control will hybridize to actin mRNA transcripts of ~2.5 kb.

Hybridization

During the hybridization phase of the protocol it is essential that the membrane is not allowed to dry completely. This will lead to irreversible binding of the probe to the membrane, high background, and will prevent successful stripping and reuse of the northern blot.

Stripping the Membrane

MessageMap northern blot membranes can be stripped and reprobed without loss of signal when performed according to the protocol (See *Stripping the Blot for Reuse*).

PROTOCOL

Preparation of Labeled Probe

Stratagene suggests radiolabeling 25–50 ng of sample DNA or RNA. The specific activity of the labeled probe should be $>1 \times 10^8$ cpm/ μ g when using fresh (less than 1 week old) radiolabel. Stratagene offers the Prime-It® II random primer labeling kit for generating high-specific activity probes. Unincorporated radiolabel must be removed from the probe prior to hybridization. Excess radionucleotides are efficiently removed from labeled probes using Stratagene's NucTrap® probe purification columns.

If using other commercially available labeling or purification kits, please follow the manufacturer's instructions. For additional labeling and purification protocols please consult references 1 and 2.

Hybridization of Radiolabeled Probe to the MessageMap® Blot

1. Generally, 180 μ l of the MiracleHyb hybridization solution per cm² of the MessageMap blot should be used for standard experiments. More or less MiracleHyb hybridization solution may be used as desired, but it is necessary to ensure that there is sufficient solution to cover the membrane at all times during the prehybridization and hybridization. The table below gives recommended minimum volumes of MiracleHyb hybridization solution for common hybridization containers.

Container*	Minimum volume of MiracleHyb™ hybridization solution
50-ml conical tubes	5 ml
heat-sealable bags or roller bottles	10 ml

2. For **double-stranded probes**, prehybridize the MessageMap blot in MiracleHyb hybridization solution at 68°C for 15 minutes.

For **oligonucleotide probes and riboprobes**, calculate the melting temperature (T_m) (see *Appendix I: Hybridization and Melting Temperatures* for the mathematical formula). Prehybridize the MessageMap blot in MiracleHyb solution for 15 minutes at 5–10°C below the T_m .

Note *It is normal for the MiracleHyb hybridization solution to become slightly opaque during the pre-hybridization step, especially at higher incubation temperatures ($\geq 50^\circ\text{C}$).*

* If 50-ml screw cap conical tubes or bottles are used, ensure that the blot is positioned with the Stratagene logo facing inward and not pressed against the walls of the vessel. If processing more than one blot, place each blot in a separate vessel

3. Aliquot into a screw-cap microcentrifuge the correct amount of MiracleHyb probe preparation buffer (50 μ l of MiracleHyb probe preparation buffer per 5 ml of prehybridization solution used in Step 2.)

Note *Always use screw-cap microcentrifuge tubes when boiling radioactive solutions*

4. Add radiolabeled probe to the MiracleHyb probe preparation buffer and mix by pipetting. For double-stranded probes, boil the diluted probe for 2 minutes. Briefly spin in a bench top microcentrifuge. **The boiling step is not necessary for oligonucleotide or RNA probes.**

For best results, use random-primed radioactive probes with the following concentration and specific activity:

Suggested Probe Concentration

1.0 $\times 10^6$ total counts/ml of hybridization solution

Specific Activity of the Probe

10⁸ cpm/ μ g or greater

5. Add the probe/MiracleHyb probe preparation buffer mixture to the prehybridization solution containing the MessageMap blot. **It is important that the probe is added to the prehybridization solution and is not pipetted directly onto the membrane.**

Hybridization Conditions

Carry out the hybridization at 68°C. (When using oligonucleotide probes or other short probes, perform the hybridization at 5–10°C below the T_m .) Probe may be hybridized to the blot (on roller bars or a similar device for gentle agitation) for as few as one to two hours for high-abundance messages or standard Southern blot analysis. For low-abundance messages or when maximum sensitivity is desired, hybridize overnight (≥ 16 hours).

Washing the Blot

Perform the following membrane washes (gentle agitation is required) for double-stranded probes, oligonucleotide probes and riboprobes:

1. Wash twice for 15 minutes each at room temperature with excess 2 \times SSC buffer and 0.1% (w/v) SDS wash solution[§].
2. Wash once for 30 minutes at 60°C with excess 0.1 \times SSC buffer and 0.1% (w/v) SDS wash solution[§] for a high-stringency wash.

[§] See *Preparation of Media and Reagents*.

Detection

Using forceps, grasp the MessageMap blot by one corner and lift it out of the wash solution. Hold it vertically so that the excess liquid is allowed to drain off. Remove the last traces of excess liquid by touching the bottom corner of the membrane to a clean paper towel or tissue. Wrap the blot in plastic wrap. Expose the wrapped MessageMap blot to autoradiography film with an intensifying screen at -80°C for 2 hours–overnight. For very low-abundance messages or increased signal intensity, film exposure may be extended up to one week without greatly increasing background.

Stripping the Blot for Reuse

1. Heat the $0.1\times$ SSC buffer and 0.1% (w/v) SDS wash solution to boiling.
2. In a glass dish, pour the $0.1\times$ SSC buffer and 0.1% (w/v) SDS wash solution over the MessageMap blot and wash the blot twice for 15 minutes each.

Proceed with the prehybridization step for the next hybridization. If the MessageMap blot will not be used immediately for a second round of hybridization, remove the excess liquid from the membrane by draining (see *Detection* above), then store the MessageMap blot in plastic wrap, desiccated.

β -actin cDNA Control

Use 25 ng of the β -actin cDNA control when preparing the control probe using random primer labeling protocols. Probes suitable for this experiment should have a specific activity of $>1 \times 10^8$ cpm/ μg . The hybridization solution for the control experiment should contain 1×10^6 cpm/ml β -actin probe. A sharp band ~ 2.5 kb* should be detected on the film after 1–2 hours at -80°C with an intensifying screen.

* In some cases the β -actin control might also hybridize to α -actin, causing two bands to appear on the blot.

APPENDIX I: HYBRIDIZATION AND MELTING TEMPERATURES

Hybridization Temperature

The hybridization temperature for oligonucleotide probes equals 5–10°C below the T_m .

Calculation of the Melting Temperature

Note *The first equation below overestimates the T_m of hybrids involving longer oligonucleotides. The second formula works only for monovalent cation concentrations of ≤ 1 M. The molar concentration of monovalent cations in MiracleHyb solution is 0.5 M.*

Oligonucleotides Shorter than 18 Bases

$$T_m = 2^\circ\text{C}(A + T) + 4^\circ\text{C}(G + C)$$

Oligonucleotides 14 Bases and Longer (up to 60–70 Nucleotides)

$$T_m = 81.5 + 16.6(\log_{10}[MC]) + 0.41(\%G + C) - (600/N)$$

where $[MC]$ is the molar concentration of monovalent cations and N is the chain length.

More discussion regarding calculation of melting temperature can be found in references 1 and 2.

TROUBLESHOOTING

Observation	Suggestion
Low hybridization signal	Check the specific activity of the labeled probe to ensure that it is $>1 \times 10^8$ cpm/ μg
	Ensure that the concentration of the probe in solution is $>1 \times 10^6$ cpm/ml of MiracleHyb hybridization solution
	Increase exposure time of blot to X-ray film
High background	Remove unincorporated nucleotides from the probe prior to hybridization
	Repeat the wash sequences
	Reduce exposure time
	Ensure that the correct amount of DNA was used in the labeling reaction

PREPARATION OF MEDIA AND REAGENTS

<p>20× SSC Stock Solution 175.3 g of NaCl (3 M final concentration) 88.2 g of sodium citrate-trisodium salt (300 mM final concentration) Adjust the pH to 7.0 with HCl Add distilled water (dH₂O) to 1 liter</p>	<p>2× SSC Buffer and 0.1% (w/v) SDS Wash Solution 100 ml of 20× SSC buffer 5 ml of 20% (w/v) SDS dH₂O to 1 liter</p>
<p>20% (w/v) SDS Stock Solution Dissolve 20 g of SDS in 90 ml of distilled water (dH₂O) Mix well and heat to 68°C if necessary Add dH₂O to a final volume of 100 ml</p>	<p>0.1× SSC Buffer and 0.1% (w/v) SDS Wash Solution 5 ml of 20× SSC buffer 5 ml of 20% (w/v) SDS dH₂O to 1 liter</p>

REFERENCES

1. Ausubel, F. M., Brent, R., Kingston, R. E., Moore, D. D., Seidman, J. G. *et al.* (1987). *Current Protocols in Molecular Biology*. John Wiley & Sons, New York.
2. Sambrook, J., Fritsch, E. F. and Maniatis, T. (1989). *Molecular Cloning: A Laboratory Manual*. Cold Spring Harbor Laboratory Press, Cold Spring Harbor, NY.

ENDNOTES

MessageMap[®], NucTrap[®] and Prime-It[®] are registered trademarks of Stratagene in the United States.

MiracleHyb is a trademark of Stratagene.